



Meaning of the different terms used in brochures

Meaning of the terms "Maximum", "Typical", "Calibrated" and "Minimum" used in EIZO Display Technologies product brochures

QUESTION: What do the terms "Maximum", "Typical", "Calibrated" and "Minimum" refer to as used in EIZO Display Technologies product brochures?

ANSWER: The specifications of medical displays are largely determined by the LCD module and backlight used. Comparing such values can be particularly difficult when different measuring techniques and assumptions are used. EIZO Display Technologies interprets some common specification terms as follows.

- **Maximum** specified value is at this level or lower (may be limited to the value at the time of delivery)
- **Typical** specified value represents a potential value that may or may not be maintained by an individual product (typical values often represent individual components and may not take external effects into account)
- **Calibrated** specified value of the display is calibrated, too (higher or lower values may be specified individually)
- **Minimum** specified value is at this level or higher (may be limited to the value at the time of delivery)